

# Glycosynthase (Endo-M-N175Q)



An enzyme that adds whole sugar chains without breaking down products.

**Increases efficiency of the sugar metastasis reaction**

**Applicable to Biosimilars**


**Making of a new functional sugar complex**




New tool for Sugar synthesis

Glycosynthase (Endo-M-N175Q) is an enzyme developed by Yamamoto, Umekawa, *et al.* through site-directed mutation of areas around the active center of Endo-M<sup>1)</sup> which is already marketed. Since the feature of Glycosynthase is efficient transglycosylation activity by using oxazoline derivatives as glycosyl donors while suppressing sugar hydrolysis activity, the resulting glycosylated products are obtained in high yield with less digestion of the products by the enzyme. Due to this feature Glycosynthase is expected to be applied as a useful tool in glyco-technology.

Fig. 1 Comparison experiment <sup>2)</sup>




Specialized in transfer activity, and inhibited hydrolysis activity



Sugar chains of clear structure

Glycoprotein authentic sample



Sugar metastasis  
Glycosynthase (Endo-M-N175Q)


The acquisition of standardized sugar complexes

Biomedicines  
Biosimilars




Umekawa and her colleagues caused transglycosylation reactions at the GlcNAc site of sperm antigen CD52 using oxazoline derivatives of the high-mannose type sugar chains or the complex type sugar chains as glycosyl donors<sup>2)</sup>. They succeeded in obtaining glycosylated products in high yield of 84% and 76%, respectively. Moreover, they also achieved transglycosylation reactions using two biologically active blood-pressure-lowering peptides, PAMP12 and Substance P, as glycosyl acceptors and an oxazoline derivative of a complex type sugar chain containing sialic acids as a glycosyl donor in 95% and 98% yield, respectively<sup>3)</sup>.

Fig. 2 Experiment example of the transglycosylation




### HPLC monitoring of the sugar metastasis reaction



Rate of metastasis  
**95%**

HPLC : ODS  
CH<sub>3</sub>CN / H<sub>2</sub>O  
UV260nm



Realizing efficient transglycosylation reactions would also be useful for expansion into glycoprotein synthesis, such as the area of biosimilars. Creation of new functional sugar complexes can be expected also.

- |              |  |                    |
|--------------|--|--------------------|
| <b>G0365</b> | <b>Glycosynthase (Endo-M-N175Q)</b><br>from <i>Mucor hiemalis</i> expressed in <i>Escherichia coli</i> (100m units/vial) | 1vial              |
| <b>A1651</b> | <b>Endo-β-N-Acetylglucosaminidase (Endo-M)</b><br>from <i>Mucor hiemalis</i> expressed in <i>Candida boidinii</i>        | 1vial              |
| <b>S0523</b> | <b>Sialylglycopeptide (SGP)</b>  | 10mg               |
| <b>D4065</b> | <b>Disialyloctasaccharide</b>  | Please contact us. |

1) K. Yamamoto, S. Kadowaki, J. Watanabe, H. Kumagai, *Biochem. Biophys. Res. Commun.* **1994**, *203*, 244.  
2) M. Umekawa, T. Higashiyama, T. Tanaka, M. Noguchi, A. Kobayashi, S. Shoda, W. Huang, L-X. Wang, H. Ashida, K. Yamamoto, *Biochim. Biophys. Acta, Gen. Subj.* **2010**, *1800*, 1203.  
3) M. Umekawa, C. Li, T. Higashiyama, W. Huang, H. Ashida, K. Yamamoto, L-X. Wang, *J. Biol. Chem.* **2010**, *285*, 511.

#### Ordering and Customer Service

#### TCI AMERICA

Tel : 800-423-8616 / 503-283-1681  
Fax : 888-520-1075 / 503-283-1987  
E-mail : Sales-US@TCIchemicals.com

#### TCI EUROPE N.V.

Tel : +32 (0)3 735 07 00  
Fax : +32 (0)3 735 07 01  
E-mail : Sales-EU@TCIchemicals.com

#### TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00  
Fax : +49 (0)6196 64053-01  
E-mail : Sales-DE@TCIchemicals.com

#### Tokyo Chemical Industry UK Ltd.

Tel : +44 (0)1865 784560  
Fax : +44 (0)1865 784561  
E-mail : Sales-UK@TCIchemicals.com

#### TCI Chemicals (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909  
Fax : 044-2262 8902  
E-mail : Sales-IN@TCIchemicals.com

#### 梯希爱(上海)化成工业发展有限公司

Tel : 800-988-0390 / 021-67121386  
Fax : 021-6712-1385  
E-mail : Sales-CN@TCIchemicals.com

#### TOKYO CHEMICAL INDUSTRY CO., LTD.

Tel : +81 (0)3-5640-8878  
Fax : +81 (0)3-5640-8902  
E-mail : globalbusiness@TCIchemicals.com

Availability, price or specification of the listed products are subject to change without prior notice. Reproduction forbidden without the prior written consent of Tokyo Chemical Industry Co., Ltd.