



乙烯化试剂


乙烯基硅烷衍生物




Vinyltrimethylsilane
25mL / 100mL
[V0067]




Triethylvinylsilane
5mL / 25mL
[T1703]




TASF: Tris(diethylamino)sulfonium Difluorotrimethylsilicate
Y. 98%
Y. Hatanaka, T. Hiyama, *J. Org. Chem.* **1988**, 53, 918.




Vinyltrimethoxysilane
25mL / 100mL / 500mL
[V0042]




Triethoxyvinylsilane
25mL / 100mL / 500mL
[V0044]




R = COMe, CO₂H, OMe, NH₂
Y. ~ 98%
A. Gordillo, M. A. Ortuño, C. L.-Mardomingo, A. Lledos, G. Ujaque, E. de Jesus, *J. Am. Chem. Soc.* **2013**, 135, 13749.




2,4,6,8-Tetramethyl-2,4,6,8-tetravinylcyclotetrasiloxane
5g / 25g
[T2523]



R = 2-NO₂, Y.85%
4-MeO, Y.86%
S. E. Denmark, C. R. Butler, *Org. Lett.* **2006**, 8, 63.




(1-Fluorovinyl)methyldiphenylsilane
1g / 5g
[F0518]




Y. 86%
T. Hanamoto, T. Kobayashi, *J. Org. Chem.* **2003**, 68, 6354.

乙烯基硼烷衍生物




2,4,6-Trivinylboroxin - Pyridine Complex
1g / 5g
[T2498]




Y. 73%
F. Kerins and D. F. O'Shea, *J. Org. Chem.* **2002**, 67, 4968.



Potassium Vinyltrifluoroborate
1g / 5g / 25g
[P1479]





R = CO₂Et or NO₂
R = CO₂Et, Y. 84%
NO₂, Y. 88%
S. Darses, G. Michaud, J.-P. Genet, *Eur. J. Org. Chem.* **1999**, 1875.




4,4,5,5-Tetramethyl-2-vinyl-1,3,2-dioxaborolane
(stabilized with Phenothiazine)

1g / 5g / 25g
[T2297]




Dibutyl Vinylboronate
(stabilized with Phenothiazine)

1g / 5g
[D3596]




其它




Tributylvinyltin

1g / 5g / 25g
[T1794]




Vinylmagnesium Bromide
(14% in Tetrahydrofuran, ca. 1mol/L)

100g
[V0053]



Vinyl Bromide
(ca. 12% in Tetrahydrofuran, ca. 1.0mol/L)

100mL
[V0127]



更多信息, 请查看我们的主页: www.TCIchemicals.com ▶▶▶

乙烯化反应

